## INgene® q PRRS UNIVERSAL

Cat. No. R.11.PRU.K.5TX/Q

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*qPCR kit for the detection of porcine reproductive and respiratory syndrome virus (PRRSV) in biological samples.*

*Real Time PCR (TaqMan™ technology). 100 reactions.*

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**Index of contents (ENGLISH)**

1.	INTRODUCTION.....	9
2.	INTENDED USE .....	9
3.	PRINCIPLE OF THE ASSAY .....	9
4.	KIT COMPOSITION .....	9
5.	CONSERVATION OF THE KIT COMPONENTS .....	9
6.	MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED.....	9
7.	WORKING GUIDELINES AND PRECAUTIONS TO AVOID CONTAMINATION .....	10
8.	SAMPLES COLLECTION AND TRANSPORT .....	10
9.	RNA EXTRACTION.....	10
10.	AMPLIFICATION OF THE GENETIC MATERIAL .....	11
11.	ANALYSIS AND INTERPRETATION OF RESULTS .....	12
12.	TROUBLESHOOTING.....	13
13.	QUALITY CONTROL.....	14
14.	TECHNICAL ASSISTANCE AND ORDERS .....	14

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9175.INGE

## 1. INTRODUCTION

**INgene q PRRS UNIVERSAL** is a real-time PCR kit suited for the detection of RNA of the porcine reproductive and respiratory syndrome virus (PRRSV) in biological samples in a simple way, with high sensitivity and specificity.

PRRSV is the causal agent of a porcine syndrome that causes reproductive failure in breeding sows and respiratory illness in young swine, causing serious productive losses. PRRSV is an enveloped virus with an RNA genome belonging to the genus *Arterivirus*, family *Arteriviridae*, order *Nidovirales*. There are two related but antigenically and genetically distinguishable strains, North American and European, with a genomic divergence of 40%. The high mutation and evolution rate of PRRSV strains causes the need to frequently update the diagnostic molecular assays to include recent genetic variants.

## 2. INTENDED USE

**INgene q PRRS UNIVERSAL** is a real-time PCR kit suited for the detection of RNA of porcine reproductive and respiratory syndrome virus (PRRSV) in biological samples of clinical interest, such as serum, lung, tonsils, oral fluids and semen.

*For research use only.*

## 3. PRINCIPLE OF THE ASSAY

**INgene q PRRS UNIVERSAL** is an assay for amplification and detection of the PRRSV genetic material based on the polymerase chain reaction (PCR). The assay consists of a real-time PCR in a one-well format. The kit contains all necessary reagents and enzymes to amplify and simultaneously detect PRRSV genetic material in extracted RNA samples. The reagents are included in two mixes that, once combined in the indicated proportions, are in optimal concentrations to detect the genetic material of PRRSV, simply by adding the extracted ARN of the problem sample. This kit contains specific primers and a TaqMan™ probe labelled with FAM fluorochrome for the specific detection of PRRSV. The primer and probe mix also includes specific primers and TaqMan™ probe labelled with Cy5 fluorochrome to amplify the Internal Positive Control (IPC) that serves as a workflow validation control. In addition, the master mix contains a passive reference fluorochrome, ROX, to allow signal normalisation in the real-time PCR platforms where this is required. The amount of ROX in the master mix is adjusted for the thermocyclers that are indicated to use LOW ROX (such as Applied Biosystems™ Quantstudio™ 5). In the case that higher ROX levels are required it should be supplemented by the end user.

The kit includes a PRRSV positive control for detection.

The kit has been tested in the Quantstudio™ 5 instrument from Applied Biosystems™ and in LightCycler® 480 II from ROCHE It has been validated in FAST mode in the instrument Quantstudio™ 5. It is compatible with other thermocyclers, if they have the appropriate fluorescence channels, but Ct values may vary among them.

## 4. KIT COMPOSITION

COMPONENTS (100 reactions)	Nº VIALS	VOLUME/VIAL
Mixture A: PRRSV and IPC specific primers and probes	2	600 µl
Mixture B: enzyme mix	2	300 µl
Positive Control A1: PRRSV amplification control	1	100 µl
Internal positive control (IPC)	1	600 µl

## 5. CONSERVATION OF THE KIT COMPONENTS

All components must be stored at -20°C immediately after arrival and until the moment of use and protected from light. The components of the kit are stable for 1 year from the manufacturing date (see expiry date on packaging).

## 6. MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Powder-free disposable gloves
- Microcentrifuge

- Tube shaker
- Real-time thermocycler
- Micropipettes (0,5-1000 µl)
- Sterile pipette tips (with filter)
- DNase/RNase-free water, molecular grade

## 7. WORKING GUIDELINES AND PRECAUTIONS TO AVOID CONTAMINATION

The following working guidelines should be carefully observed:

- Use DNase/RNase-free disposable consumables.
- Use DNase/RNase-free water, molecular grade.
- Use sterile filtered tips.
- Ensure good homogenisation of the kit components once defrosted.
- During use, keep all the kit components on ice/cooling block.
- Avoid repeated cycles of freezing and defrosting of the kit components, as it can reduce kit performance. Protect kit reagents from light exposure until their immediate use.

Special precautions to avoid ambient or cross-contaminations that could cause false positive results:

- The positive PCR control should be stored separately from the rest of kit reagents.
- Store the kit components in dedicated areas separated from samples or amplification products.
- Use separate working areas with dedicated material and unidirectional workflow, for every step of the process: (1) sample manipulation and extraction, (2) qPCR setup, (3) qPCR amplification and (4) just in the case needed, analysis of amplified products.
- If the above recommendation is not feasible, separate at least two separated areas, pre-PCR (steps 1 and 2) and post-PCR (steps 3 and 4).
- During the qPCR setup, the positive control should be manipulated and dispensed in the last place, avoiding the contact with the other kit components.

## 8. SAMPLES COLLECTION AND TRANSPORT

**INgene q PRRS Universal** is a real-time PCR kit suited for the detection of RNA of porcine reproductive and respiratory syndrome virus (PRRSV) in biological samples of clinical interest, such as serum, lung, tonsils, oral fluids and semen.

Samples should be collected and handled carefully to avoid contamination and should arrive in the testing laboratory as soon as possible, following national regulations for their transport. The samples must be cooled from the moment they are obtained and during the transport until they are processed. Processing needs to occur within 24 hours after obtaining the samples. In case a longer time is envisaged, freezing of samples is recommended.

## 9. RNA EXTRACTION

Any extraction procedure which yields good quality of the extracted material could be used. Several commercial suppliers could provide high quality nucleic-acid extraction, following the manufacturer instructions, such as: Applied Biosystems™, ROCHE, QIAGEN®, MACHEREY-NAGEL™, others.

Recommended controls during extraction step:

**Internal Positive Control:** The kit includes an Internal Positive Control (IPC) to be added in the lysis buffer, in amount of 5 µl per sample. It should not be added directly in the samples. This IPC is based on a heterologous sequence that is purified and finally amplified in every PCR reaction, simultaneously to the pathogen nucleic acid. The IPC follows the same workflow as the samples and allows to detect false negatives.

**Extraction Control:** Parallel to the extraction of the genetic material of the samples, it is advised to perform a mock-extraction, in which the biological sample is replaced by nuclease-free water, to check possible contamination during the extraction procedure.

**CAUTION:** *If amplification is not going to be done immediately, preserve the RNA at -20°C or -80°C.*

**10. AMPLIFICATION OF THE GENETIC MATERIAL**

REQUIRED MATERIALS

- Ice bucket or cooling block
- RNA extracted from the samples
- Mixture A and B – KEEP ON ICE AND PROTECTED FROM LIGHT meanwhile PCR setup, store at -20°C afterwards
- Positive amplification control A1 (PRRSV)
- DNAase-RNAase-free water (molecular grade)
- Tubes/plates with optical caps/sealing compatible with thermocycler
- DNAase-RNAase-free tubes (1,5 ml)

PROCEDURE

1. Prepare as many tubes or PCR plate wells for the amplification as samples to be processed, adding additional tubes for the positive and the negative amplification controls. Identify positions avoiding direct label on the tubes/plates at it could interfere with optical reading.
2. Thaw and keep mixes A and B on ice/cooling block. Make sure that they are correctly homogenised before pipetting. Prepare a fresh amplification mix by adding in one tube the amount needed for the total number of samples to be processed, including the necessary controls. Keep the mixes A and B (and the amplification mix) always on ice/cooling block and protected from light.
3. Prepare the amplification mix on ice as follows. It is advisable to prepare an excess amount of mix (calculate an extra 10% for all reagents) in order to compensate for possible losses of volume during pipetting:

	Mixture A	Mixture B	Master Mix
Per sample	10 µl	5 µl	15 µl
For 10 samples	100 µl	50 µl	150 µl

4. Once mixture is prepared, homogenise it correctly. Place the PCR plate or strips on ice/cooling block and add 15 µl of the final reaction mix to each well.
5. Add 5 µl of each extracted sample to the corresponding tube or reaction well of the PCR plate, or alternatively, of each of the reaction controls indicated.

Type of sample/control	Component	Volume
Sample	Extracted DNA from biological samples	5 µl
Extraction control	Extraction of a water control, in parallel to the samples	5 µl
Positive control (C+)	Positive amplification control A1 (PRRSV)	5 µl
Negative control (C-)	Nuclease-free water	5 µl

6. Close/seal the tubes/plate and mix carefully the content, make sure the liquid is at the bottom of the tubes/wells. Otherwise, spin down briefly in a centrifuge.
7. Insert the plate/tubes in the real-time instrument and follow the specific instructions to programme the reaction according to the following conditions:

	Temperature (°C)	Time	Cycles
Reverse Transcription	50	10 min	1
Denaturalization	95	2 min	1
Amplification	95	5 sec	40
	60 *	30 sec	

**FAST mode:** the kit has been also validated in FAST mode with the cycler Quantstudio™ 5 (Applied Biosystems™).

8. Reading of fluorescence takes place during the hybridisation/elongation step (asterisk in the table\*), using the indicated fluorescence channels:

	Reporter	Quencher
PRRSV	FAM	None
IPC	Cy5	None
Passive reference	ROX	

## 11. ANALYSIS AND INTERPRETATION OF RESULTS

At the end of the run, results should be analysed with the software of the thermocycler, following the manufacturer instructions. From each sample, fluorescence data originated from channel FAM is obtained. It is recommended in first place to select the automatic analysis mode for threshold and baseline determination. Alternatively, in the case that the automatic threshold is set incorrectly, adjust it manually following the thermocycler instructions. The threshold should be placed above the background and in the initial phase of the exponential curve. The baseline starts and end cycles, should be set before any significant fluorescence is detected.

Interpretation of fluorescence data:

FAM channel: pathogen detection. A positive result should show a typical amplification curve with a Ct value < 40. A negative result implies absence of amplification (undetected) or a Ct ≥ 40.

Cy5 channel: Internal Positive Control detection. A positive result should show a typical amplification curve with a Ct value <40. Commonly registered Ct values are between 30 < Ct < 40.

In all channels: A negative result implies absence of amplification (undetected) or Ct ≥ 40.

### Criteria for the kit validation:

**The assay will be considered as valid when the C+ have a Ct value within the range 25 ± 4 and the negative control Ct ≥ 40 in the FAM channel.**

Before interpreting the samples result, check the validation criteria of the controls:

Control	FAM channel Pathogen detection	Cy5 channel Internal Positive Control	Results
Negative control (C-)	undetected (ND), Ct ≥ 40	Not relevant	Valid PCR
	Ct < 40	Not relevant	Invalid PCR (see section 12.1)
Positive control (C+)	Ct 25 ± 4	Not relevant	Valid PCR
	Ct < 21 o Ct > 29	Not relevant	Invalid PCR (see section 12.2)
Extraction control	undetected (ND), Ct ≥ 40	Ct < 40	Valid extraction control
	undetected (ND), Ct ≥ 40	Ct ≥ 40 (ND)	Invalid extraction (see section 12.3)
	Ct < 40	Not relevant	Invalid extraction (see section 12.4)

In the case that the validation criteria of the controls are met, the interpretation of the sample results should be as follows:

**Table:** Interpretation of sample results.

Case	Sample	Detection channels		Result	Interpretation
		FAM	Cy5		
1	Problem sample	Positive (Ct < 40)	Positive (Ct < 40) or Negative (Ct ≥ 40)	Positive	Sample contains RNA of PRRSV
2	Problem sample	Negative (Ct ≥ 40)	Positive (Ct < 40)	Negative	RNA of PRRSV is NOT detected
3	Problem sample	Negative (Ct ≥ 40)	Negative (Ct ≥ 40)	Invalid	Possible inhibition (see section 12.5)

## 12. TROUBLESHOOTING

12.1. Detection of fluorescence signal on FAM channels (pathogen) in the negative PCR control.

Possible cause	Solution
Contamination during PCR setup	<ul style="list-style-type: none"> <li>Repeat PCR reaction with new reagents.</li> <li>If tubes are being used, close each one immediately after adding the sample.</li> <li>As a strict rule, the positive control must be manipulated and added in last place and avoiding any contact with the rest of reaction components.</li> <li>Ensure that the workplace, equipment and all instruments used for performing the PCR are decontaminated, i.e. free of nucleic acids.</li> </ul>

12.2. Absence of fluorescence signal on FAM and/or HEX channel in the positive PCR controls.

Possible cause	Solution
Wrong fluorochrome channel selected	<ul style="list-style-type: none"> <li>Check that FAM channel is selected for all samples and controls to be analysed.</li> </ul>
Wrong analysis settings (automatic or manually selected)	<ul style="list-style-type: none"> <li>Verify the adequate positioning of the threshold and baseline, and the corrected selection of passive fluorescence dye (ROX) in the case needed, according to the instrument requirements.</li> </ul>

12.3. Absence of fluorescence signal on Cy5 channel in the negative extraction control.

Possible cause	Solution
Extraction efficiency can differ among extraction kits. The CIP used has been tested with: Maxwell® RSC blood DNA kit (PROMEGA) and MagMax Pathogen RNA/DNA kit (Thermo Fisher Scientific)	<ul style="list-style-type: none"> <li>Title the CIP for the kit in use.</li> <li>Use one of the kits with which the CIP has been validated.</li> </ul>

12.4. Detection of fluorescence signal on FAM channel in the negative extraction control.

Possible cause	Solution
Contamination occurred during the extraction process	<ul style="list-style-type: none"> <li>Repeat the RNA extraction process and PCR using new reagents.</li> <li>Ensure that the workplace, equipment and all instruments used for performing PCR are decontaminated (free of nucleic acids).</li> </ul>

12.5. Absence of fluorescence signal in the three fluorescence channels FAM, HEX and Cy5 in the case of negative samples.

Possible cause	Solution
PCR inhibition	<ul style="list-style-type: none"> <li>Dilute 1/40 the extracted RNA material and repeat the PCR reaction.</li> <li>In the case that inhibition persists with the diluted sample, repeat the RNA extraction from the biological sample.</li> </ul>

**13. QUALITY CONTROL**

Each production batch of this kit has been validated under our ISO-9001 and 14001 certified Quality Management System.

**14. TECHNICAL ASSISTANCE AND ORDERS**

For further information concerning the assay and its performance:

**GOLD STANDARD DIAGNOSTICS MADRID, S.A**  
 e-mail: [biologiamolecular.Spain@eu.goldstandarddiagnostics.com](mailto:biologiamolecular.Spain@eu.goldstandarddiagnostics.com)

For placing orders, please contact:

e-mail: [cs.spain@eu.goldstandarddiagnostics.com](mailto:cs.spain@eu.goldstandarddiagnostics.com)  
**GOLD STANDARD DIAGNOSTICS MADRID, S.A.**  
 Calle de los Hermanos García Noblejas, 39-8º  
 28037, Madrid, España.  
 Tel: +34 91 368 0501

Websites: [www.goldstandarddiagnostics.es](http://www.goldstandarddiagnostics.es) [www.goldstandarddiagnostics.com](http://www.goldstandarddiagnostics.com)